

Rodent Anesthesia System



Features

- Maintain accurate flow rate to any circuit, unlike stopcock or gang-valve delivery systems
- Simple to use - no calculations
- Comes with short, easy to follow protocol
- Specific flow to each animal, regardless of number of circuits in use or length of tubing
- Connect breathing circuits or chambers - your choice
- With one vaporizer and 6 individual connections, connect as many circuits or chambers as you need in any combination
- Lower preset flows mean less waste gases created, resulting in lower trace exposure to investigator
- Preset flows at the flip of a switch - No flowmeter adjustments
- Tabletop, Mobile, or Wall-mounted

RC2 - Rodent Circuit Controller System includes:

vaporizer; 2 - induction chambers; 4 - breathing circuits; 4 - nose cones; 22' - 22mm evacuation tubing; and 6 charcoal filters.

Info: Fabrizio Barbieri

KF Technology Srl
(+39) 06 454.34.179 - info mobile:(+39) 339 533.03.22 - fax.(+39) 06 9725.3131- SkyPe : kft2002
fabrizio@kftechnology.it - www.kftechnology.it

Masks & Replacement Diaphragms

These high-quality masks feature a standard 15mm male connector, clear cone for full visualization and a highly flexible, replaceable rubber diaphragm.



Info: Fabrizio Barbieri

KF Technology Srl
(+39) 06 454.34.179 - info mobile:(+39) 339 533.03.22 - fax.(+39) 06 9725.3131- SkyPe : kft2002
fabrizio@kftechnology.it - www.kftechnology.it