

Programmable Peristaltic Pump NE-9000



Description:

- Have a need to accurately fill multiple bottles, vials or other containers with a liquid? The learn and repeat mode of the NE-9000 Peristaltic Pump from New Era Pump systems makes even the smallest of businesses highly productive right from day one & MORE

Features:

- Dispense and withdrawal pumping
- Dispense and withdrawal volumes separately accumulated
- Programmable dispense volumes
- Non-volatile memory of all operating parameters and pumping program
- Programmable Phases allowing complex pumping applications and interaction with external devices.
- RS-232 bi-directional control
- Worldwide power supplies available
- Built-in pump network driver. Pump network supports up to 100 pumps and other devices
- Two modes of RS-232 control, Basic and Safe. Safe mode provides communication error detection, loss of communication detection, and automatic transmitting of alarm conditions.
- TTL I/O with software filtered control inputs to eliminate glitches and ringing on the control inputs.
- Configurable TTL operational trigger.
- Power Failure Mode: Restarts the Pumping Program after a power interruption
- Audible Alarm