

Tweezers & Forceps

DUMONT #5

500085

Material: Dumostar - non-magnetic, non-corrosive
Length: 11cm Tips: 0.025 x 0.015 mm
US \$ 68



500233

Material: Dumostar - non-magnetic, non corrosive
Length: 11 cm Tips: 0.10 x 0.06 mm
US \$ 44



14098

Material: Dumoxel - non-magnetic, non-corrosive
Length: 11cm Tips: 0.1 x 0.06 mm
US \$ 31



500342

Material: Inox
Length: 11 cm Tips: 0.10x0.06 mm
US \$ 25



500341

Material: Inox, Biology
Length: 11 cm Tips: 0.05x0.01 mm
US \$ 30



14095

Material: Stainless Steel, Medical Biology
tweezers
Length: 11cm Tips: 0.05 x 0.01 mm
US \$ 36



501201

Material: Stainless Steel, Medical, 90° bent
Length: 11 cm Tips: 0.05 x 0.01 mm
US \$ 36



	Hardness (Rockwell)	Max. Temp. Resistance	Corrosion Resistance	Magnetic
Stainless Steel & Inox	55-56	350°C	Excellent	Yes
Dumoxel	36	350°C	Excellent	No
Dumostar	62	550°C	100% Non-Corrosive	No
Titanium	37	550°C	100% Non-Corrosive	No

14101

Material: Stainless Steel, Medical Biology, 45°
Length: 11cm Tips: 0.05 x 0.01 mm
US \$ 39



500232

Material: Stainless Steel, Medical Biology
Length: 11 cm Tips: curved, 0.02 x 0.06 mm
US \$ 40



500232

Material: Stainless Steel, Medical Biology
Length: 11 cm Tips: curved, 0.02 x 0.06 mm
US \$ 40



14096

Material: Titanium
Length: 11cm Tips: 0.1 x 0.06 mm
US \$ 38



DUMONT #7

14188

Material: Dumoxel - non-magnetic, non-corrosive
Length: 12 cm Tips: 0.17 x 0.1 mm
\$ 39



501302

Material: Stainless Steel
Length: 11cm Tips: 0.17 x 0.10mm serrated jaws
\$ 56



14097

Material: Inox
Length: 12 cm long Tips: 0.17 x 0.1 mm
\$ 32



14187

Material: Titanium - non-magnetic, non-corrosive
Length: 12 cm Tips: 0.17 x 0.1 mm
\$ 40



Dressing Forceps

500092

Adson Forceps
Length: 12 cm, 1x2 teeth
\$ 12

Tip Profile:



14226
Adson Forceps
Length: 12 cm, serrated tips
\$ 17



14142
Graefe Forceps
Length: 7 cm, serrated straight 0.7 mm tips
\$ 21

Tip Profile:



14140
Graefe Forceps
Length: 7 cm, 1x2 teeth, straight 0.7 mm tips
\$ 21

Tip Profile:



500363
Dressing Forceps
Length: 15.5 cm, serrated (pictured)
\$ 18

500364
Dressing Forceps
Length: 25.5 cm, serrated
\$ 18

500365
Dressing Forceps
Length: 30.5 cm, serrated
\$ 18



501217
Dressing Forceps
Length: 12.5 cm, serrated
\$ 18

15914
Iris Forceps
Length: 10 cm, serrated straight 0.8 mm tips
\$ 17

Tip Profile:



15916
Iris Forceps
Length: 10 cm, 1x2 teeth, straight 0.8 mm tips
\$ 17

Tip Profile:



15917
Iris Forceps
Length: 10 cm, 1x2 teeth, curved 0.8 mm tips
\$ 17

Tip Profile:



500455

Filter Forceps (great for gel handling)
Polypropylene finger grip, length: 11 cm, straight flat jaw
\$ 22




500451

Micro Hemostatic Forceps
Length: 12.5 cm long, serrated, 0.4 x 0.8 mm tips
\$ 39

500452

Micro Hemostatic Forceps
Length: 12 cm long, serrated, 0.4 x 0.8 mm tips
\$ 39



501240

Mixer Hemostatic Forceps "Baby", right angle
Length: 14 cm
\$ 33



Rochester-Oschner Hemostatic Forceps

501709 -- Straight, length 14 cm \$ 22

501710 -- Curved, length 14 cm \$ 22

501326 -- Straight, length 16 cm \$ 23

501711 -- Curved, length 16 cm \$ 23



Scissors

14177

McPherson-Vannas Scissors
Length: 7 cm, straight 3 mm blades, 0.1 mm tips
\$ 146



501778

Vannas Scissors, Super Fine
Length: 8 cm, straight 3 mm blades, 0.015 X 0.015 mm tips
\$ 194

